

4ti-0955
96 Well Semi-Skirted PCR Plate, Roche Style

Brooks Life Sciences are committed to providing innovative products that advance applications and workflows for our customers.

We strive to continuously improve processes and efficiency to ensure our customers receive the service they expect from Brooks Life Sciences.

Unfortunately, due to an unexpected issue with tooling, we are currently experiencing a small variation from our accepted standards for 4ti-0955. This variation does not affect the form or function of 4ti-0955 but we have observed a less than adequate fit with our Cap Strips 4ti-0751 and 4ti-0783.

Customers using any other seals, such as 4ti-0560, will be unaffected by this variation. If your application requires the use of Cap Strips, please contact your territory manager who will be able to suggest alternative products.

We apologize for the impact that this has on your work flow and can assure you that Brooks Life Sciences is assessing all options to improve this.

If you have any questions regarding the status of your order, please contact your territory manager or email:

4ti-Orders@brooks.com



Tracy Andrews | Senior Product Manager | Consumables and Instruments | 4titude

Phone: +44 (0) 1306 884 885 | Mobile: +44 (0) 7966 307379

Email: tracy.andrews@brooks.com | Web: www.brookslifesciences.com