

Silicone Sealing Cap Mat Compatibility Table

Silicone Sealing Cap Mats	Microplates											
	96 Round Well Microplate, Low Binding 200 µl, V-shaped bases (4ti-LB0109)	96 Round Well Microplate 300 µl, U-shaped bases (4ti-0110)	96 Round Well Microplate 330 µl, U-shaped bases (4ti-0116)	96 Round Well Microplate 330 µl, V-shaped bases (4ti-0117)	96 Round Deep Well Storage Microplate 1.2 ml, U-shaped bases (4ti-0120)	96 Round Deep Well Storage Microplate 2.0 ml, U-shaped bases (4ti-0130)	96 Round Deep Well Storage Microplate, For Magnetic Separators 1 ml, V-shaped bases (4ti-0125)	96 Round Deep Well Storage Microplate, For Magnetic Separators, Low Binding 1 ml, V-shaped bases (4ti-LB0125)	96 Square Deep Well Storage Microplate 1.2 ml, U-shaped bases (4ti-0126)	96 Square Deep Well Storage Microplate 2.2 ml, V-shaped bases (4ti-0132)	96 Square Deep Well Storage Microplate 2.2 ml, U-shaped bases (4ti-0136)	384 Square Deep Well Storage Microplate 190 µl, V-shaped bases (4ti-0147)
96 Round Well, clear, pierceable (4ti-0124)							✓					
96 Round Well, white, non-pierceable (4ti-0135)					✓							
96 Square Well, clear, pierceable (4ti-0137)								✓	✓	✓		
96 Round Well, clear, pierceable (4ti-0138)	✓		✓	✓		✓	✓	✓				
384 Square Well, clear, pierceable (4ti-0139)												✓