






NOTE:
ALL FACES TO HAVE A SMOOTH GLOSS FINISH (HIGHLY POLISHED ON TOOL)
UNDIMENSIONED FEATURES BY 3D MODEL



ACTUAL SIZE

ISSUE	DATE	AMENDMENT



WWW.4ti.CO.UK

	CUSTOMER		DRAWN	ZOLTAN E.		TITLE CrystalCaps (8-strip caps) MARKETING DRAWING	
	CUSTOMER DRG. No.		ISSUE	DATE	12.03.2015		
DO NOT SCALE	MATERIAL POLYCARBONATE - Panlite L 1225Z 100		CHECKED	PAUL D.			
IF IN DOUBT ASK	TOLERANCES		DATE		DRAWING No. 4ti-0755		
	GENERAL TOLERANCES (UNLESS OTHERWISE STATED) 0 ±0.20 0.0 ±0.10 0.00 ±0.05 ANGULAR ±0°30' (ALL DIMENSIONS IN MILLIMETRES)		SCALE 5:1		SHEET No. 1 OF 1		ISSUE A